

SC-3008-DL Sequential Circulator Pneumatic Compression Device

The ultimate gradient sequential compression device for the treatment of Lymphedema, Venous Insufficiency, and Wound Healing. The SC-3008-DL Sequential Circulator is as unique as the patient population it services. All patients have different needs so our custom developed technology measures and then delivers pressure to each and every chamber in the garment to give the best possible treatment. The SC-3008-DL lets the user program the pressure setting in any and all of the 8 chambers, always keeping greater pressure distal so as to avoid any reflux. Complete with a compliance meter to monitor patient use, the SC-3008-DL also has a countdown clock and a timer that shuts off the pump after one hour.




Bio Compression Systems, Inc.
120 West Commercial Avenue
Moonachie, NJ 07074

Phone: 201-939-0716
Fax: 201-939-4503
E-mail: biosystems@biocompression.com
www.biocompression.com

Quick Start Guide

SC-3008-DL





Key Function

1. POWER ON/OFF BUTTON
2. MODE BUTTON
3. INCREASE BUTTON
4. DECREASE BUTTON
5. LED DISPLAY
6. CHAMBER NUMBER LED
7. RECEPTOR PORTS FOR TUBING LATCH CONNECTORS

Setting the Pressure

- Turn on pump by pressing the Power button.
- You have 5 seconds after turning on the pump to press and hold the Up and Down arrows.
- After 5 seconds of holding the arrows, the LED will light up with 50 and 1. This is the pressure and chamber number.



- Let go of arrow buttons and adjust chamber 1 pressure by using the arrow buttons.
- Press the "Mode" button to move to chamber 2 and adjust pressure using the arrow buttons.
- Continue Step 5 to adjust pressure in all chambers.
- Press the Power button to lock in settings and begin treatment.

Setting the Mode

- Turn on pump by pressing the Power button.
- You have 5 seconds after turning on the pump to press and hold the "Mode" button.

- After 5 seconds of holding the "Mode" button, "1 Hr" will appear in the window indicating that you are in the countdown timer mode. The pump is shipped in this mode already, if you are happy with a 1 hour countdown timer you do not have to set the mode.



- If you would like to use the countdown timer mode but for a different amount of time, use the Up and Down arrows to adjust the timer between 10 and 120 minutes. Press the Power button to lock in the time, this will be set for all future treatments.

- If you would like to run the pump in Continuous mode, after step 3 press the "Mode" button. You will now see "CO" in the LED screen. Press the Power button and your pump will always run in Continuous mode.



- If you would like to run the pump with Pre-Therapy, after Step 3 press the "Mode" button 2 times. The LED display will now show "PrE 0". The "0" indicates that Pre-Therapy is currently turned off. Press the Up arrow one time and you will now see "Pre 1". The "1" indicates that Pre-Therapy is on. Hit the "Mode" button 1 time to go back to the countdown timer mode, 2 times to go to the continuous mode, then hit the Power button to lock the settings in place. Note- you can not lock the settings in place while "Pre 0" or "Pre 1" are showing.



Hard Reset to return pump to factory Setting

- The pump must be off, do not hit the Power button.
- Press and hold the Up arrow, Down arrow, and Mode button all at the same time for 5 seconds.
- The LED will light up after 5 seconds, you can now release the buttons.
- Press the "Mode" button one time and the LED will turn off.



Reminders

- Your pump will always arrive in the 1 Hour Countdown mode and pressure set to 50 mmHg with a 1 mmHg gradient sequential setting.
- While running the pump in the countdown timer mode, you can press and hold the "Mode" button to toggle between pressure settings and remaining time.