

120 West Commercial Avenue, Moonachie, New Jersey 07074  
Phone (800) 888-0908 / Alt. Phone (201) 939-0716 / Fax (201) 939-4503  
www.biocompression.com

## USING YOUR SEQUENTIAL CIRCULATOR IN CHINA, KOREA, AND WESTERN EUROPE!

### VOLTAGE ADAPTER INFORMATION

"I have successfully used the Model #3004 and #3008 Sequential Gradient Circulator in China, Korea, and Western Europe utilizing the Radio Shack Model 1/85 85 Watt Step-Down Converter. It will convert 220/240 VAC to 110/120 VAC. It and the circulator work fine with 50 Hz. This converter is not off the shelf at most Radio shacks, they have to order it. Also, the user should get the appropriate electrical outlet adapters depending on where they travel. All necessary information is in the catalog in Radio Shack's retail stores.

I would also recommend taking a small roll of duct tape or other strong tape. After connecting the converter, with or without an electrical connection adapter, to the wall outlet, it may become physically unstable due to it's weight and need taping to the wall.

The adapter is small, so traveling with it is no problem. It is approximately 2.5in. x 2.5in. x 4.125in."

*-As written by a satisfied user*



**\*Other sources for travel supply are:** MAGELLAN's @ 800-962-4943 & INTER POWER @ 800-662-2290.