

BIOCRYO GARMENT SIZING CHART

Upper Extremity

	<u>Axilla</u>	<u>Bicep</u>	<u>Elbow</u>	<u>Forearm</u>	<u>Wrist</u>	<u>Axilla to Fingertip Length</u>
Arm GC-3035-L	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5	35.5"/91.2
Arm w/ Shoulder GC-3035-SH	27"/69.6	26"/66.0	24"/61.0	23"/58.4	19.5"/49.5	31+8"/78.7+20.3

Lower Extremity

	<u>Thigh below Crotch</u>	<u>Mid Thigh</u>	<u>Mid Knee</u>	<u>Mid Calf</u>	<u>Ankle</u>	<u>Arch (Instep)</u>	<u>Crotch to Heel</u>
Half Leg GC-3045-H			29"/73.7	25.5"/64.8	24.5"/62.2	21.5"/54.6	
Small Leg GC-3045-S	32"/81.3	31"/78.7	25"/63.5	24.5"/62.2	24.5"/62.2	21.5"/54.6	27"/69.6
Medium Leg GC-3045-M	32"/81.3	31"/78.7	27"/68.6	25"/63.5	24.5"/62.2	21.5"/54.6	32"/81.3
Large Leg GC-3045-L	32"/81.3	31"/78.7	29"/73.7	25.5"/64.8	24.5"/62.2	21.5"/54.6	36"/91.4

