






BIO VEST APPLICATION

Using Velcro type attachment system, after pulling over head, snugly pull front together and attach.



BIO VEST APPLICATION WITH BOTH UNILATERAL & BILATERAL UPPER EXTREMITY SLEEVES

Using Velcro type attachment circumferentially around top of sleeve, simply attach to arm portion of vest.



120 West Commercial Avenue, Moonachie, NJ 07074

Phone: 201-939-0716 / Toll-Free: 800-888-0908 / Fax: 201-939-4503
E-mail: biosystems@biocompression.com / www.biocompression.com

