

# Bio Compression Systems

## *Elite 8 BioVest Sizing Chart*

		Circumferences					Lengths			
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips	Thorax	
GV-3010-S-L	Small Left	63.5	59.6	53.3	50.8	48.2	64.7	74.9	Adjustable to 140	71
GV-3010-S-R	Small Right	63.5	59.6	53.3	50.8	48.2	64.7	74.9	Adjustable to 140	71
GV-3010-M-L	Medium Left	63.5	59.6	53.3	50.8	48.2	71.1	81.2	Adjustable to 140	71
GV-3010-M-R	Medium Right	63.5	59.6	53.3	50.8	48.2	71.1	81.2	Adjustable to 140	71
GV-3010-L-L	Large Left	63.5	59.6	53.3	50.8	48.2	78.7	88.9	Adjustable to 140	71
GV-3010-L-R	Large Right	63.5	59.6	53.3	50.8	48.2	78.7	88.9	Adjustable to 140	71

**U.S. BASED  
MANUFACTURER**

*NOTE: Measurements Above are in Centimeters and are the maximum allotted point on limb*