

Bio Compression Systems

Elite 8 BioVest Sizing Chart

		Circumferences					Lengths		
Part Number	Size	Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips	Thorax
GV-3010-S-L	Small Left	25"	23.5"	21"	20"	19"	25.5"	29.5"	Adjustable 28" to 55"
GV-3010-S-R	Small Right	25"	23.5"	21"	20"	19"	25.5"	29.5"	Adjustable 28" to 55"
GV-3010-M-L	Medium Left	25"	23.5"	21"	20"	19"	28"	32"	Adjustable 28" to 55"
GV-3010-M-R	Medium Right	25"	23.5"	21"	20"	19"	28"	32"	Adjustable 28" to 55"
GV-3010-L-L	Large Left	25"	23.5"	21"	20"	19"	31"	35"	Adjustable 28" to 55"
GV-3010-L-R	Large Right	25"	23.5"	21"	20"	19"	31"	35"	Adjustable 28" to 55"

**U.S. BASED
MANUFACTURER**

NOTE: Measurements Above are in **Inches** and are the maximum allotted point on limb