

Bio Compression Systems

Lower Extremity/ 4-Chamber Garment Sizing Chart

Part Number	Size	Circumferences						Lengths		
		Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease	
GS-3045-H	Half	Standard	NA		73.7	68.6	59.7	55.9	48.3	NA
GW-3045-H		Wide			77.5	71.1	73.7	58.4		
GA-3045-H		Adjustable			81.3	76.2	59.7	55.9		
GWA-3045-H		Adj. Wide			85.1	78.7	73.7	58.4		
GN-3045-S	Small	Narrow	64.8	63.5	55.9	50.8	45.7	49.5	64.8	68.6
GS-3045-S		Standard	78.7	76.2	71.1	62.2	54.6	49.5		
GW-3045-S		Wide	90.2	86.4	77.5	71.1	64.8	54.6		
GA-3045-S		Adjustable	86.4	83.8	78.7	69.9	54.6	49.5		
GWA-3045-S		Adj. Wide	97.8	94	85.1	78.7	64.8	54.6		
GN-3045-M	Medium	Narrow	69.9	64.8	58.4	50.8	45.7	49.5	71.1	81.3
GS-3045-M		Standard	81.3	76.2	71.1	62.2	54.6	49.5		
GW-3045-M		Wide	92.7	86.4	76.2	71.1	64.8	57.2		
GXW-3045		Extra Wide	111.8	104.1	91.4	86.4	76.2	57.2	73.7	81.3
GA-3045-M		Adjustable	88.9	86.4	78.7	69.9	54.6	49.5	71.1	81.3
GWA-3045-M		Adj. Wide	100.3	92.7	85.1	78.7	64.8	57.2		
GXWA-3045		Adj. Extra Wide	119.4	111.8	99.1	94	76.2	57.2	73.7	81.3
GN-3045-L	Large	Narrow	71.1	66	61	53.3	45.7	49.5	81.3	91.4
GS-3045-L		Standard	83.8	81.3	76.2	66	54.6	49.5		
GW-3045-L		Wide	94	91.4	81.3	73.7	61	52.1		
GA-3045-L		Adjustable	94	91.4	83.8	73.7	54.6	49.5		
GWA-3045-L		Adj. Wide	101.6	96.5	88.9	81.3	61	57.2		

U.S. BASED
MANUFACTURER

NOTE: Measurements above are in Centimeter and are the maximum allotted point on limb