

Bio Compression Systems

Lower Extremity/ 4-Chamber Garment Sizing Chart

Bio Compression Systems										
Lower Extremity/ 4-Chamber Garment Sizing Chart										
Circumferences										
Lengths										
Part Number	Size	Thigh below Crotch	Mid Thigh	Mid Knee	Mid Calf	Ankle	Arch (Instep)	Crotch to Heel (Length)	Inguinal Crease	
GS-3045-H	Half	Standard	NA		29"	27"	23.5"	22"	19"	NA
GW-3045-H		Wide			30.5"	28"	29"	23"		
GA-3045-H		Adjustable			32"	30"	23.5"	22"		
GWA-3045-H		Adj. Wide			33.5"	31"	29"	23"		
GN-3045-S	Small	Narrow	25.5"	25"	22"	20"	18"	19.5"	25.5"	27"
GS-3045-S		Standard	31"	30"	28"	24.5"	21.5"	19.5"		
GW-3045-S		Wide	35.5"	34"	30.5"	28"	25.5"	21.5"		
GA-3045-S		Adjustable	34"	33"	31"	27.5"	21.5"	19.5"		
GWA-3045-S		Adj. Wide	38.5"	37"	33.5"	31"	25.5"	21.5"		
GN-3045-M	Medium	Narrow	27.5"	25.5"	23"	20"	18"	19.5"	28"	32"
GS-3045-M		Standard	32"	30"	28"	24.5"	21.5"	19.5"		
GW-3045-M		Wide	36.5"	34"	30"	28"	25.5"	22.5"		
GXW-3045		Extra Wide	44"	41"	36"	34"	30"	22.5"	29"	32"
GA-3045-M		Adjustable	35"	34"	31"	27.5"	21.5"	19.5"	28"	32"
GWA-3045-M		Adj. Wide	39.5"	36.5"	33.5"	31"	25.5"	22.5"		
GXWA-3045		Adj. Extra Wide	47"	44"	39"	37"	30"	22.5"	29"	32"
GN-3045-L	Large	Narrow	28"	26"	24"	21"	18"	19.5"	32"	36"
GS-3045-L		Standard	33"	32"	30"	26"	21.5"	19.5"		
GW-3045-L		Wide	37"	36"	32"	29"	24"	20.5"		
GA-3045-L		Adjustable	37"	36"	33"	29"	21.5"	19.5"		
GWA-3045-L		Adj. Wide	40"	38"	35"	32"	24"	22.5"		

U.S. BASED
MANUFACTURER

NOTE: Measurements above are in Inches and are the maximum allotted point on limb