

## Bio Compression Systems

### Upper Extremity/ 4-Chamber Garment Sizing Chart

Part Number	Size	Circumferences					Lengths	
		Axilla	Bicep	Elbow	Forearm	Wrist	Axilla to Fingertips	Shoulder to Fingertips
GS-3035-S	Small	59.7	55.9	49.5	44.5	38.1	60.9	67.3
GS-3035-M	Medium	68.6	66.0	61.0	53.3	45.7	71.1	78.7
GS-3035-L	Large	64.8	61.0	54.6	48.3	40.6	78.7	90.2
GS-3035-SH/S	Small	59.7	55.9	49.5	44.5	38.1	60.9	67.3 +20.3
GS-3035-SH/M	Medium	68.6	66.0	61.0	53.3	45.7	71.1	78.7 +20.3
GS-3035-SH/L	Large	64.8	61.0	54.6	48.3	40.6	78.7	90.2 +20.3

**U.S. BASED  
MANUFACTURER**

*NOTE: Measurements Above are in Centimeters and are the maximum allotted point on limb*