

Bio Compression Systems

BioPants Size Chart

Inches	Circumferences							Length					Part Number
	Upper Abdomen	Lower Abdomen	Groin	Mid thigh	Knee	Mid Calf	Ankle	Foot/ Instep	Abdomen (GT)	Full Leg (AG)	Lower Leg (AD)	Foot	
Small	44"	44"	30"	29"	28"	25"	21"	20"	15"	29"	18"	12"	GBA-3045-S-2
Medium	44"	44"	31"	30"	29"	26"	22"	20"	16"	32"	20"	12"	GBA-3045-M-2
Large	48"	48"	35"	34"	32"	29"	24"	20"	18"	35"	22"	12"	GBA-3045-L-2

NOTE: Measurements above are in Inches and are the maximum allotted point on limb

**U.S. BASED
MANUFACTURER**