

120 West Commercial Avenue, Moonachie, NJ 07074  
Phone (201) 939-0716 / Fax (201) 939-4503  
www.biocompression.com

## USING YOUR SEQUENTIAL CIRCULATOR IN CHINA, KOREA, AND WESTERN EUROPE!

### VOLTAGE ADAPTER INFORMATION

"I have successfully used the Model #3004 and #3008 Sequential Gradient Circulator in China, Korea, and Western Europe utilizing the Enercell Model V85 Watt Step-Down Converter. It will convert 220/240 VAC to 110/120 VAC. It and the circulator work fine with 50 Hz. This converter can be found on the internet for purchase. Also, the user should get the appropriate electrical outlet adapters depending on where they travel.

The adapter is small, so traveling with it is no problem. It is approximately 2.5in. x 2.5in. x 4.125in."



-As written by a satisfied user

*Please Note: The Model SC-3004-DL & the Model SC-3008-DL CANNOT be used with a converter. Please call our office for more information. Thank you!*

**\*Other sources for travel supply are:** MAGELLAN's @ 800-962-4943 & INTER POWER @ 800-662-2290.